

# MedSmart Programming and Event Reporting Worksheet



Complete this worksheet and call the **Product Support and Programming Center at 1-800-853-1458** to have your MedSmart activated and programmed. Please verify that all information supplied for device programming is complete, accurate and current.

<b>MedSmart Serial # (located on bottom of dispensing unit):</b> _ _ _ _ _											
<b>MedSmart User Data</b>											
Name				DOB				Gender			
Address											
City				State				Zip			
Home Phone			Cell Phone			Email					
<b>Device Programming Information</b>											
Time Zone (circle):    Eastern    Central    Mountain    Pacific    Alaska    Hawaii											
# Doses/Alarms per Day (circle):			1	2	3	4	5	6			
Time of Alarms:	1		2		3		4		5		6
		am pm		am pm		am pm		am pm		am pm	
<b>Primary Caregiver Data</b>											
Name				Relation							
Address											
City				State				Zip			
Home Phone			Cell Phone			Email					
<b>Call User On Phone After Missed Dose Alert is Sent</b>			<input type="checkbox"/> Yes		<input type="checkbox"/> No		If yes, select which phone to call				
			<input type="checkbox"/> Home Phone		<input type="checkbox"/> Cell Phone						
<b>Alert/Message Rules and Delivery Time</b>											
<b>Alert/Message</b>			<b>Number of Occurrences Prompting an Alert</b>								
<b>Missed Dose Alert</b>			1	2	3	4	5	6			
			<i>How many doses are missed before an alert is sent?</i>								
<b>No Device Call-In Alert</b>			1	2	3	4	5	6			
			<i># of days in row where device does not call-in before alert is sent?</i>								
<b>Tray Low/Empty Message</b>			1	2	3	4	5	6			
			<i>How many days of pills remain in tray before alert is sent?</i>								
<b>Next Day Alert/Message Delivery Time:</b> ___:___ am (delivered in user's time zone)											
<b>Product Support and Programming Center at 1-800-853-1458</b>											

Alert Recipient Data		
<b>Recipient 1 - <input type="checkbox"/> (check if same as caregiver) - If yes just fill out alert information</b>		
Name		Relation
Home Phone	Cell Phone	Email
Alert Method (check one): <input type="checkbox"/> Text Message (cell phone) <input type="checkbox"/> Email <input type="checkbox"/> Voice Message		
<b>Recipient 2 -</b>		
Name		Relation
Home Phone	Cell Phone	Email
Alert Method (check one): <input type="checkbox"/> Text Message (cell phone) <input type="checkbox"/> Email <input type="checkbox"/> Voice Message		
<b>Recipient 3</b>		
Name		Relation
Home Phone	Cell Phone	Email
Alert Method (check one): <input type="checkbox"/> Text Message (cell phone) <input type="checkbox"/> Email <input type="checkbox"/> Voice Message		
<b>Recipient 4</b>		
Name		Relation
Home Phone	Cell Phone	Email
Alert Method (check one): <input type="checkbox"/> Text Message (cell phone) <input type="checkbox"/> Email <input type="checkbox"/> Voice Message		

**IMPORTANT WARNING**

The MedSmart unit is not an emergency response system. The device is intended solely to serve as a medication reminder and dispensing system and is not in any way intended to be relied upon to provide any healthcare or medical assistance or to insure that the pills will be properly dispensed and taken by the user.